

Anti-HNRNPd antibody (1-306) (STJ118125)
STJ118125

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/IP/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 IP:0.5 Mu g-4 Mu g antibody for 200 Mu g-400 Mu g extracts of whole cells ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specif
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	3184
Gene Symbol	HNRNPd
Uniprot ID	HNRNPd_HUMAN
Immunogen	
Immunogen Region	1-306
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-306 of human HNRNPd (NP_002129.2).
Immunogen Sequence	MSEEQFGGDGAAAAATAAVG GSAGEQEGAMVAATQGAAAA AGSGAGTGGGTASGGTEGGS AESEGAKIDASKNEEDEGHs NSSPRHSEATAQREEWKMF IGGLSWDTTKDLKDYFSKF GEVVDCTLKLPITGRSRGF GFVLKESSESVDKVMQDQKEH KLNGKVIDPKRAKAMKTKEP VKKIFVGGGLSPDTPEEKIRE YFGGFGEVESIELPMDNKTN KRRGFCFTFKEEEPVKKI



Immunohistochemistry of paraffin-embedded rat liver using HNRNPd antibody (STJ118125) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human lung cancer using HNRNPd antibody (STJ118125) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using HNRNPd antibody (STJ118125) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081