

Anti-SERPINA4 antibody (178-427) (STJ118148)

STJ118148

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5267
Gene Symbol	SERPINA4
Uniprot ID	KAIN_HUMAN
Immunogen	
Immunogen Region	178-427
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 178-427 of human Kallistatin (SERPINA4) (NP_006206.2).
Immunogen Sequence	TIQLNDHVKKETRGKIVDL VSELKKDVLMLVNIYFKA LWEKPFISSRTTPKDFYVDE NTTVRVPMMLQDQEHHWYHL DRYLPCSVLRMDYKGDATVF FILPNQGKMREIEEVLTPEM LMRWNNLLRKRNFYKKLELH LPKFSISGSYVLDQILPRLG FTDLFSKWADLSGITKQKQL EASKSFHKATLDVDEAGTEA AAATSFAIKFFSAQTNRHIL RFNRPFLVIFSTSTQSVL