

Anti-FSTL1 antibody (1-308) (STJ118234)

STJ118234

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	11167
Gene Symbol	FSTL1
Uniprot ID	FSTL1_HUMAN
Immunogen	
Immunogen Region	1-308
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-308 of human FSTL1 (NP_009016.1).
Immunogen Sequence	MWKRWLALALVAVAWVRA EEELRSKSKICANVFCGAGR ECAVTEKGPTCLCIEQCKP HKRPVCGSNGKTYLNHCELH RDACLTGSKIQVDYDGHCKE KKSVPSPASPVCYQSNRDE LRRRIQWLEAEIIPDGWFS KGSNYSEILDKYFKNFDNGD SRLDSSSEFLKFVEQNETAIN ITTYPDQENKLLRLGLCVDA LIELSDENADWKLSTFQELK CLNPSFNPPEKKCALEDET




Immunofluorescence analysis of Human placenta using FSTL1 Rabbit polyclonal antibody (STJ118234) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of NIH/3T3 cells using FSTL1 Rabbit polyclonal antibody (STJ118234) at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of C6 cells using FSTL1 Rabbit polyclonal antibody (STJ118234) at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using FSTL1 Rabbit polyclonal antibody (STJ118234) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081