

Anti-COPG2 antibody (542-871) (STJ118262)
STJ118262

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:100 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	26958
Gene Symbol	COPG2
Uniprot ID	COPG2_HUMAN
Immunogen	
Region	542-871
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 542-871 of human COPG2 (NP_036265.3).
Immunogen Sequence	LNATYIFNGLTVSVPGMEKA LHQYTLPESEKPFDMKSIPL AMAPVFEQKAEITLVATKPE KLAPSRQDIFQEQLAAIPEF LNIGPLFKSSEPVLTEAET EYFVRCIKHMFNTNHIVQFD CTNTLNDQLLEKVTVMQEPS DSYEVLSCIPAPSLPYNQPG ICYTLVRLPDDPTAVAGSF SCTMKFTVRDCDPNTGVPDE DGYDDEYLEDLEVTVSDHI QKVLKPNFAAAWEEVGDTF



Immunofluorescence analysis of U-2 OS cells using COPG2 Rabbit polyclonal antibody (STJ118262) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of COPG2 in paraffin-embedded mouse brain using COPG2 Rabbit polyclonal antibody (STJ118262) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using COPG2 Rabbit polyclonal antibody (STJ118262) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJ5000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081