

Anti-PIK3R4 antibody (1-330) (STJ118273)

STJ118273

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	IHC-P:1:50-1:100 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	30849
Gene Symbol	PIK3R4
Uniprot ID	P13R4_HUMAN
Immunogen	
Immunogen Region	1-330
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-330 of human PIK3R4/VPS15 (NP_055417.1).
Immunogen Sequence	MGNQLAGIAPSQILSVESYF SDIHDFEYDKSLGSTRFFKV ARAKHREGLVVKVFAIQDP TLPLTSYKQELEELKIRLNS AQNCLPFQKASEKASEKAAM LFRQYVRDNLVDRISTRPFL NNIEKRWIAFQILTAVDQAH KSGVRHGDIKTENVMVTSWN WVLLTDFASFKPTYPEDNP ADFNYFFDTSRRRTCYIAPE RFVDGGMFATELEYMRDPST PLVDLNSNQRTRGELKRAM



Immunofluorescence analysis of NIH/3T3 cells using PIK3R4/VPS15 Rabbit polyclonal antibody (STJ118273) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using PIK3R4/VPS15 antibody (STJ118273) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 180s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081