

Anti-NUDCD1 antibody (304-583) (STJ118364)

STJ118364

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84955
Gene Symbol	NUDCD1
Uniprot ID	NUDC1_HUMAN
Immunogen	
Immunogen Region	304-583
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 304-583 of human NUDCD1 (NP_116258.2).
Immunogen Sequence	IQFLPDHINIVLKDHFLEG KLYSSIDHESSTWIKESNS LEISLIKKNEGLTWPELVIG DKQGEIIRDSAQCAIAERL MHLTSEELNPNPDKEKPPCN AQEECDIFFEESSSLCRF DGNTLKTTHWNLGNSNYLF SVIVDPKEMPCFCLRHVDVA LLWQPHSSKQDDMWEHIATF NALGYVQASKRDKKFFACAP NYSYAALCECLRRVFIYRQP APMSTVLYNRKEGRQVGQV



Western blot analysis of various lysates using NUDCD1 Rabbit polyclonal antibody (STJ118364) at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081