

Anti-ZNF766 antibody (60-190) (STJ118370)

STJ118370

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90321
Gene Symbol	ZNF766
Uniprot ID	ZN766_HUMAN
Immunogen	
Immunogen Region	60-190
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 60-190 of human ZNF766 (NP_001010851.1).
Immunogen Sequence	MMKQRTEPWTVENEMKVAKN PDRWEGIKDINTGRSCAVRS KAGNKPITNQLGLTFQLPLP ELEIFQGEKGIYECNQVQKF ISHSSSVSPLQRIYSGVKTH IFNKHNRNDFVDFPLLSQEQK AHIRRKPYECN



Western blot analysis of various lysates using ZNF766 Rabbit polyclonal antibody (STJ118370) at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081