

Anti-PLEKHA1 antibody (285-404) (STJ118449)

STJ118449

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	59338
Gene Symbol	PLEKHA1
Uniprot ID	PKHA1_HUMAN
Immunogen	
Immunogen Region	285-404
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 285-404 of human PLEKHA1 (NP_067635.2).
Immunogen Sequence	SGAIVAQRGPGRSASSEHPP GPSESKHAFRPTNAATATSH STASRSNSLVSTFTMEKRGF YESLAKVKPGNFKVQTVSPR EPASKVTEQALLRPQSKNGP QEKDCDLVLDLDDASLPVSDV