

Anti-C1QTNF7 antibody (24-190) (STJ118471)

STJ118471

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	114905
Gene Symbol	C1QTNF7
Uniprot ID	C1QT7_HUMAN
Immunogen	
Immunogen Region	24-190
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-190 of human C1QTNF7 (NP_001128642.1).
Immunogen Sequence	QPRGNQLKGENYSPRYICSI PGLPGPPGPPGANGSPGPHG RIGLPGRDGRDGRKGEKGEK GTAGLRGKTGPLGLAGEKGD QGETGKKGPIGPEGEKGEVG PIGPPGPKGDRGEQGDPLP GVCRCGSIVLKSAFSVGITT SYPEERLPPIFNKVLFNAGE HYNPATG