

Anti-GP1BA antibody (19-259) (STJ118488)
STJ118488

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2811
Gene Symbol	GP1BA
Uniprot ID	GP1BA_HUMAN
Immunogen	
Immunogen Region	19-259
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-259 of human CD42b (NP_000164.5).
Immunogen Sequence	ICEVSKVASHLEVNC DKRNL TALPPDLPKD TTILH LSEN L LYTFSLATLMPYTRLTQLNL DRCELTKLQVDGTLPLVGLTLD SHNQLQSLPLLGGQTLPAL TVLDVSNRLTSLPLGALRG LGELQELYLKGNELKTLPPG LLTPTPKLEKLSLANNNLTE LPAGLLNGLENLDTLLQEN SLYTIPKGFFGSHLLPFAFL HGNPWLCNCEILYFRWLQD NAENVYVWKQGV DVKAMTS



Immunofluorescence analysis of THP-1 cells using GP1BA Polyclonal Antibody (STJ118488) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081