

Anti-SLC18A3 antibody (473-532) (STJ118508)

STJ118508

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6572
Gene Symbol	SLC18A3
Uniprot ID	VACHT_HUMAN
Immunogen	
Immunogen Region	473-532
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 473-532 of human SLC18A3 (NP_003046.2).
Immunogen Sequence	GLLTRSRSEKDVLLDEPPQG LYDAVRLRERPVSGQDGEPR SPPGPFDAEDDYNYYTTS



Immunofluorescence analysis of paraffin-embedded mouse brain using SLC18A3 Rabbit polyclonal antibody (STJ118508) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of SLC18A3 in paraffin-embedded rat brain using SLC18A3 Rabbit polyclonal antibody (STJ118508) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.