

Anti-USO1 antibody (663-962) (STJ118519)
STJ118519

GENERAL INFORMATION





Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8615
Gene Symbol	USO1
Uniprot ID	USO1_HUMAN
Immunogen	
Immunogen Region	663-962
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 663-962 of human USO1 (NP_003706.2).
Immunogen Sequence	DLQLEELRQVSTLKCQNEQ LQTAVTQQVSIQHQHKDQYN LLKIQLGKDNQHQGSYSEGA QMNGIQPEEIGRLREEIEEL KRNQELLQSQLTEKDSMIEN MKSSQTSQTNEQSSAIVSAR DSEQVAELKQELATLKSQLN SQSVEITKLQTEKQELLQKT EAFKSVVEVQGETETIATK TTDVEGRLSALLQETKELKN EIKALSEERTAIKEQLDSSN STIAILQTEKDKLELEITD



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081