

Anti-PROZ antibody (220-340) (STJ118522)

STJ118522

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8858
Gene Symbol	PROZ
Uniprot ID	PROZ_HUMAN
Immunogen	
Immunogen Region	220-340
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 220-340 of human PROZ (NP_003882.1). SLLHRNITVKTYFNRTSQDP LMIKITHVHVHVMRYDADAGE NDLSLLEWPIQCPGAGLP VCTPEKDFAEHLLIPRTRGL LSGWARNGTDLGNSLTTRPV TLVEGEECGQVLNVTVTTRT Y
Immunogen Sequence	



Western blot analysis of various lysates using PROZ Rabbit polyclonal antibody (STJ118522) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081