

Anti-CAB39 antibody (1-341) (STJ118565)

STJ118565

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:1000
Range	IF/CC: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51719
Gene Symbol	CAB39
Uniprot ID	CAB39_HUMAN
Immunogen	
Immunogen Region	1-341
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-341 of human CAB39 (NP_057373.1).
Immunogen Sequence	MPFPFGKSHKSPADIVKNLK ESMVLEKQDISDKKAEKAT EEVSKNLVAMKEILYGTNEK EPQTEAVAQLAQELYNSGLL STLVADLQLIDFEGKKDVAQ IFNNILRRQIGTRTPVVEYI CTQQNILFMLLKGYESPEIA LNCGIMLRECIRHEPLAKII LWSEQFYDFFRYVEMSTFDI ASDAFATFKDLLTRHKLLSA EFLEQHYDRFFSEYEKLLHS ENYVTKRQSLKLLGELLDD



Immunofluorescence analysis of NIH/3T3 cells using CAB39 rabbit polyclonal antibody (STJ118565) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using CAB39 rabbit polyclonal antibody (STJ118565) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081