

Anti-VIRMA antibody (1400-1700) (STJ118641)

STJ118641

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	IHC-P:1:50-1:200 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	25962
Gene Symbol	VIRMA
Uniprot ID	VIR_HUMAN
Immunogen	
Region	1400-1700
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1400-1700 of human KIAA1429 (NP_056311.2).
Immunogen Sequence	FMRQILNSDTIGCCGDDNGL MEVEGAHTSRTMSINAAELK QLLQSKEESPENLFLELEKL VLEHSKDDDNLDLLDSVVG LKQMLESSGDPLPLSDQDVE PVL SAPESLQNLFNRTAYV LADVMDDQLKSMWFTPFQAE EIDTDLDLVKVDLIELSEKC CSDFDLHSELSERFLSEPSS PGRITKTTKGFLGKHKHETI ITSSGKSEYIEPAKRAHVVP PPRGRGRGGFGQGIRPHDI



Immunohistochemistry analysis of paraffin-embedded mouse kidney using KIAA1429 Rabbit polyclonal antibody (STJ118641) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human esophageal cancer using KIAA1429 Rabbit polyclonal antibody (STJ118641) at dilution of 1:200 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Western blot analysis of extracts of Mouse lung, using KIAA1429 antibody (STJ118641) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081