

Anti-A2M antibody (1155-1474) (STJ118764)

STJ118764

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	IHC-P: 1:50-1:200
	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2
Gene Symbol	A2M
Uniprot ID	A2MG_HUMAN
Immunogen	
Immunogen Region	1155-1474
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1155-1474 of human Alpha-2-Macroglobulin (Alpha-2-Macroglobulin (A2M)) (NP_000005.2).
Immunogen Sequence	ALAGNQDKRKEVLKSLNEEA VKKDNSVHWERPQKPKAPVG HFYEPQAPSAEVEVMTSYVLL AYLTAPAPTSSEDLTSATNI VKWITKQQNAQGGFSSTQDT VVALHALSKYGAATFTRTGK AAQVTIQSSGTFSSKFQVDN NNRLLLQQVSLPELPGEYSM KVTGEGCVYLQTSCLKYNILP EKEEFPFALGVQTLPTQCDCE PKAHTSFQISLSVSYTGSRAS ASNMAIVDKMVSQFIPLK



Immunohistochemistry analysis of Alpha-2-Macroglobulin (A2M) in paraffin-embedded human placenta using Alpha-2-Macroglobulin (A2M) Rabbit polyclonal antibody (STJ118764) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081