

Anti-ARAF antibody (371-606) (STJ118774)

STJ118774

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:200
	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	369
Gene Symbol	ARAF
Uniprot ID	ARAF_HUMAN
Immunogen	
Immunogen Region	371-606
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 371-606 of human ARAF (NP_001645.1).
Immunogen Sequence	GFMTTRPGFAITQWCEGSSL YHHLHVADTRFDMVQLIDVA RQTAQGM DY L HAKNIIHRDL KSNIFLHEGLTVKIGDFGL ATVKTRWSGAQP LEQPSGSV LWMAAEVIRMQDPNPYSFQS DVYAYGVWLYELMTGSLPYS HIGCRDQIIFMVGRGYLSPD LSKISSNCPKAMRRLSDCL KFQREERPLFPQILATIELL QRS LPKIERSASEPSLHRTQ ADELPA CLLSAARLVP



Immunofluorescence analysis of HeLa cells using ARAF antibody (STJ118774) at dilution of 1:100. Blue: DAPI for nuclear staining.