

Anti-KRT13 antibody (269-458) (STJ118821)
STJ118821

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:100 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	3860
Gene Symbol	KRT13
Uniprot ID	K1C13_HUMAN
Immunogen	269-458
Region	
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 269-458 of human Cytokeratin 13 (Cytokeratin 13 (KRT13)) (NP_705694.2).
Immunogen Sequence	DLTRVLAEMREQYEAMAERN RRDAEEWFHAKSAELNKEVS TNTAMIQTSKTEITELRRTL QGLEIELQSQLSMKAGLENT VAETECRYALQLQIQLIS SIEAQLSELRSEMECQNEY KMLLDIKTRLEQEIATYRSL LEGQDAKMIGFPSSAGSVSP RSTSVTTTSSASVTTTSSNAS GRRTSDVRRP




Immunofluorescence analysis of human placenta using Cytokeratin 13 (Cytokeratin 13 (KRT13)) Polyclonal Antibody (STJ118821) at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.



Immunohistochemistry analysis of paraffin-embedded human stomach using Cytokeratin 13 (Cytokeratin 13 (KRT13)) antibody (STJ118821) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human breast cancer using Cytokeratin 13 (Cytokeratin 13 (KRT13)) antibody (STJ118821) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using Cytokeratin 13 (Cytokeratin 13 (KRT13)) antibody (STJ118821) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081