

Anti-NMT1 antibody (1-250) (STJ118836)

STJ118836

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IHC-P:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID	4836
Gene Symbol	NMT1
Uniprot ID	NMT1_HUMAN
Immunogen	
Immunogen Region	1-250
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human NMT1 (NP_066565.1).
Immunogen Sequence	MADESETAVKPPAPPLPQMM EGNGNGHEHCSDCENEEDNS YNRGGLSPANDTGAKKKKKK QKKKKEGSETDSAQDQPVK MNSLPAERIQEIQAIELFS VGQGPAKTMEEASKRSYQFW DTQVPKLGVEVNTHTGPVEP DKDNIRQEPYTLPGGFTWDA LDLGDRGVKLKELYTLNENY VEDDDNMFRFDYSPELLWA LRPPGWLPQWHCGVRVSSR KLVGFISAIPANIHIYDTE



Immunohistochemistry of paraffin-embedded human placenta using NMT1 antibody (STJ118836) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse kidney using NMT1 antibody (STJ118836) at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse heart using NMT1 antibody (STJ118836) at dilution of 1:100 (40x lens).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081