

Anti-BCLAF1 antibody (330-552) (STJ118901)

STJ118901

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9774
Gene Symbol	BCLAF1
Uniprot ID	BCLF1_HUMAN
Immunogen	
Immunogen Region	330-552
Specificity	A synthetic peptide corresponding to a sequence within amino acids 330-552 of human BCLAF1 (NP_001070909.1).
Immunogen Sequence	TAKTGKFLKSPPLHKNLDAR EKSTFREESPLRIKMIASDS HRPEVKLKMAPVPLDDSNRP ASLTKDRLLASTLVHSVKE QEFRSIFDHIKLPQASKSTS ESFIQHIVSLVHHVKEQYFK SAAMTLNERFTSYQKATEEH STRQKSPEIHRRIDISPSTL RKHTRLAGEERVFKENQKG DKKLRCDSADLRHDIDRRRK ERSKERGDSKGSRESSGSRK QEK



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081