

Anti-LZTS1 antibody (337-596) (STJ118923)

STJ118923

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:1000-1:5000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11178
Gene Symbol	LZTS1
Uniprot ID	LZTS1_HUMAN
Immunogen	
Immunogen Region	337-596
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 337-596 of human LZTS1 (NP_066300.1).
Immunogen Sequence	QLQQEKRLRQELSLMKEQ DLLETKLRSYEREKTSFGPA LEETQWEVCQKSGEISLLKQ QLKESQTEVNAKASEILGLK AQLKDRTRGKLEGLRLTQDL EGALRTKGLEVCENELQR KKNEAELLREKVNLLQEQLQ ELRAQAALARDMGPPTFPED VPALQRELERLRAELREERQ GHDQMSSGFQHERLVWKEEK EKVIQYQKQLQQSYVAMYQR NQRLEKALQQLARGDSAGE



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081