

Anti-TREM1 antibody (1-234) (STJ118962)

STJ118962

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	54210
Gene Symbol	TREM1
Uniprot ID	TREM1_HUMAN
Immunogen	
Immunogen Region	1-234
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human TREM1 (NP_061113.1).
Immunogen Sequence	MRKTRLWGLLWMLFVSELRA ATKLTEEKYELKEGQTLQDK CDYTLKFASSQKAWQIIRD GEMPKTLACTERPSKNSHPV QVGRILEDYHDHGLLRVRM VNLQVEDSGLYQCVIYQPPK EPHMLFDRIRLVTKGFSGT PGSNENSTQNVYKIPPTTTK ALCPLYTSPRTVTQAPPKST ADVSTPDSEINLTNVTDIIR VPFVNIULLAGGFLSKSLV FSVLFAVTLRSFVP



Western blot analysis of various lysates using TREM1 Rabbit polyclonal antibody (STJ118962) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081