

Anti-RHBG antibody (1-100) (STJ118985)

STJ118985

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	57127
Gene Symbol	RHBG
Uniprot ID	RHBG_HUMAN
Immunogen	
Immunogen Region	1-100
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RHBG (NP_065140.3).
Immunogen Sequence	MAGSPSRAAGRRLQLLLCL FLQGATAVLFAVFVRYNHKT DAALWHRSNHSNADNEFYFR YPSFQDVHAMVFGFGFLMV FLQRYGFSSVGFTFLAFA