

Anti-HRNR antibody (1-200) (STJ119047)

STJ119047

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	388697
Gene Symbol	HRNR
Uniprot ID	HORN_HUMAN
Immunogen Region	1-200
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human HRNR (NP_001009931.1).
Immunogen Sequence	MPKLLQGVITVIDVYQYAT QHGEYDTLNKAEKELLENE FHQILKNPNDPDTVDIILQS LDRDHNKKVDFTEYLLMIFK LVQARNKIIGKDYCQVSGSK LRDDTHQHQQEEQEETEKEEN KRQESSFSHSSWSAGENDSY SRNVRGSLKPGTESISRRLS FQRDFSGQHNSYSGQSSSYG EQNSDSHQSSGRGQCQSGSG