

Anti-GATA6 antibody (70-160) (STJ119058)

STJ119058

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2627
Gene Symbol	GATA6
Uniprot ID	GATA6_HUMAN
Immunogen	
Immunogen Region	70-160
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 70-160 of human GATA6 (NP_005248.2).
Immunogen Sequence	AAGPPARSLLLSSYASHPFG APHGSPAGVAGPGGNLSSW EDLLFTDLQQAATASKLLW SSRGAKLSPFAPEQPEEMYQ TLAALSSQGPA



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081