

Anti-PLA2G5 antibody (1-138) (STJ119123)

STJ119123

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5322
Gene Symbol	PLA2G5
Uniprot ID	PA2G5_HUMAN
Immunogen	
Immunogen Region	1-138
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human PLA2G5 (NP_000920.1).
Immunogen Sequence	MKGLPLAWFLACSPAVQG GLLDLKSMIEKVTGKNALTN YGFYGCYCGWGGRGTPKDGTDWCCWAHDHCYGRLEEKGCN IRTQSYKYRFAWGVVTCEPG PFCHVNLCACDRKLVYCLKR NLRSYNPQYQYFPNILCS