

Anti-ACVR1B antibody (24-126) (STJ119153)

STJ119153

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	91
Gene Symbol	ACVR1B
Uniprot ID	ACV1B_HUMAN
Immunogen	
Immunogen Region	24-126
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-126 of human ACVR1B (NP_004293.1).
Immunogen Sequence	SGPRGVQALLCACTSCLQAN YTCETDGACMVSIFNLDGME HHVRTCIPKVELVPAGKPFY CLSSEDLRNTHCCYTDYCNR IDLRVPSGHLKEPEHPSMWG PVE