

Anti-APOF antibody (220-300) (STJ119159)

STJ119159

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	319
Gene Symbol	APOF
Uniprot ID	APOF_HUMAN
Immunogen	
Immunogen Region	220-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 220-300 of human APOF (NP_001629.1).
Immunogen Sequence	DLGYDLLMTMAGMSGGPMGL AISAALKPALRSGVQQLIQY YQDQKDANISQPETTKEGLR AISDVSDLEETTTLASFISE V