

Anti-DNMT3A antibody (500-700) (STJ119195)

STJ119195

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1788
Gene Symbol	DNMT3A
Uniprot ID	DNM3A_HUMAN
Immunogen	
Immunogen Region	500-700
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human DNMT3A (NP_715640.2).
Immunogen Sequence	SVTQKHQEWGPFDLVIGGS PCNDLSIVNPARKGLYEGTG RLFFEFYRLLHDARPKEGDD RPFFWL FENVVAMGVSDKRD ISRFLSENPVIMIDAKEVSAA HRARYFWGNLPGMNRPLAST VNDKLELQCELEHGRIAKFS KVRTITTRSNSIKQGKDQHF PVFMNEKEDILWCTEMERVF GFPVHYTDVSNMSRLARQL L



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081