

Anti-TNFSF11 antibody (165-317) (STJ119274)

STJ119274

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IHC-P/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8600
Gene Symbol	TNFSF11
Uniprot ID	TNF11_HUMAN
Immunogen	
Immunogen Region	165-317
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 165-317 of human TNFSF11 (NP_003692.1).
Immunogen Sequence	FAHLTINATDIPSGSHKVSLSWYHNRGWAKISNMTFSNGKLIVNQDGFYLYANICFRH HETSGDLATEYLQLMVYVTKTSIKIPSSHTLMKGGSTKYW SGNSEFHFYSINVGGFFKLRSGEEISIEVSNPSLLDPDQD ATYFGAFKVRDID



Immunohistochemistry analysis of TNFSF11 in paraffin-embedded human placenta using TNFSF11 Rabbit polyclonal antibody (STJ119274) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of TNFSF11 in paraffin-embedded human tonsil using TNFSF11 Rabbit polyclonal antibody (STJ119274) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081