

Anti-SAP30 antibody (C-Term) (STJ119278)

STJ119278

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8819
Gene Symbol	SAP30
Uniprot ID	SAP30_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human SAP30 (NP_003855.1).
Immunogen Sequence	IELDKSARHLYICDYHKNLI QSVRNRRKRKGSDDGGDSP VQDIDTPEVDLYQLQVNTLR RYKRHFKLPTPLNKAQLV EIVGCHFRSIPVNEKDTLTY FIYSVKNDKNKSDLKVDGSGV H