

Anti-Kctd19 antibody (1-300) (STJ119418)

STJ119418

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	IF/ICC:1:50-1:200
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	279499
Gene Symbol	Kctd19
Uniprot ID	KCD19_MOUSE
Immunogen	
Immunogen Region	1-300
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of mouse Kctd19 (NP_001288102.1). MEEPGGGLHESAEDLFHFNVG GWHFSVPRSKLAQFPDSSLW KEASALTSENQRLFIDRDG STFRHVHYLYTSKLSFSSC
Immunogen Sequence	AELNLLYEQALGLQMLPLQ TLDNLKEGRHHLRVRPADIP VAERASLNYWRTWKCSKPS DFIKSPAFTGLHDKAPLGL MDTPLLDTEEEVHYCFPLD LVAKYPSLVTEDNLLWLAET VALIECECSEFRFVNFRLS HKILLPDNFSNIDVLEAEV



Immunofluorescence analysis of paraffin-embedded mouse testis using Kctd19 Rabbit polyclonal antibody (STJ119418) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.