

Anti-RNF14 antibody (125-474) (STJ119441)

STJ119441

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9604
Gene Symbol	RNF14
Uniprot ID	RNF14_HUMAN
Immunogen	
Immunogen Region	125-474
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 125-474 of human RNF14 (NP_004281.1).
Immunogen Sequence	AWMQFLKEETLAYLNIVSPF ELKIGSQKKVQRRTAQASPN TELDFGGAAGSDVDQEEIVD ERAVQDVESLSNLIQEILDF DQAQIQKCFNSKFLCSICF CEKLGSECMYFLECRHVVYCK ACLKDYFEIQIRDGQVQCLN CPEPKCPSVATPGQVKELVE AELFARYDRLLQSSLDLMA DVVYCPRPCCQLPVMQEPGC TMGICSSCNFAFCTLRLTY HGVSPCKVTAEKLMDLRNE