

Anti-ABHD6 antibody (51-336) (STJ119552)

STJ119552

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB: 1:500-1:2000
Range	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	66082
Gene Symbol	Abhd6
Uniprot ID	ABHD6_MOUSE
Immunogen	
Immunogen Region	51-336
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 51-336 of human ABHD6 (NP_079617.2).
Immunogen Sequence	RYAHHEDYQFCYSFRGRPGH KPSILMLHGFSAHKDMWLSV VKFLPKNLHLVCVDMPGHEG TTRSSLDDLSIVGQVKRIHQ FVECLKLNKKPFHLIGTSMG GHVAGVYAAYYPSDVCSLSL VCPAGLQYSTDNPFVQRLKE LEESAAIQKIPLIPSTPEEM SEMLQLCSYVRFKVPQQILQ GLVDVRIPHNSFYRKLFLFI VNEKSRYSLHENMDKIKVPT QIIWQKQDQVLDVSGADIL