

Anti-MAP1S antibody (830-1059) (STJ119569)

STJ119569

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:1000
Range	ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55201
Gene Symbol	MAP1S
Uniprot ID	MAP1S_HUMAN
Immunogen	
Immunogen Region	830-1059
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 830-1059 of human MAP1S (NP_060644.4).
Immunogen Sequence	VPPPLDPSSICMVDPEMLP PKTARQTENVSRTRKPLARP NSRAAAPKATPVAAAKTKGL AGGDRASRPLSARSEPSEKG GRAPLSRKSSSTPKTATRGPS GSASSRPGVSATPPKSPVYL DLAYLPSGSSAHLVDEEFFQ RVRALCYVISGQDQRKEEGM RAVLDAALLASKQHWDRDLQV TLIPTFDSVAMHTWYAETHA RHQALGITVLGSNSMVMQD DAFFACKVEF