

**Anti-CLCN6 antibody (1-80) (STJ119574)**  
STJ119574**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	
<b>Applications</b>	ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	Lot specific
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
<b>Formulation</b>	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	1185
<b>Gene Symbol</b>	CLCN6
<b>Uniprot ID</b>	CLCN6_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	1-80
<b>Specificity</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human CLCN6 (NP_001277.1).
<b>Immunogen Sequence</b>	MAGCRGSLCCCRWCCCGE RETRTPEELTILGETQEEED EILPRKDYESLDYDRCINDP YLEVLETMDNKKGRRYEAVK