

Anti-FAM127B antibody (1-113) (STJ119724)

STJ119724

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	26071
Gene Symbol	RTL8A
Uniprot ID	RTL8A_HUMAN
Immunogen	
Immunogen Region	1-113
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human FAM127B (NP_001071640.1).
Immunogen Sequence	MDGRVQLMKALLAGPLRPAA RRWRNPFPETFDGDTDL PEFIVQTSSYMFVDENTFSN DALKVTLITRLTGPALQWV IPYIRKESPLLNDYRGFLAE MKRVFGWEDEEDF



Immunofluorescence analysis of L929 cells using FAM127B Rabbit polyclonal antibody (STJ119724) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using FAM127B Rabbit polyclonal antibody (STJ119724) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.