

Anti-FASTKD5 antibody (495-764) (STJ119787)

STJ119787

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/IF/ICC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution	WB:1:500-1:2000
Range	IF/ICC:1:50-1:200
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60493
Gene Symbol	FASTKD5
Uniprot ID	FAKD5_HUMAN
Immunogen	
Immunogen Region	495-764
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 495-764 of human FASTKD5 (NP_068598.1).
Immunogen Sequence	SPGFVRLAQERTKFDLLKEL YTLDTGVGIECPDYRGNRLS THLQQEGSELLWYLAEKDMN SKPEFLETVFLLLETMLGGPQ YVKHHMILPHTRSSDLEVQL DVNLKPLPFNREATPAENVA KLRLEHVGVSLTDDL MNKLL KGKARGHFQGGKTESEPGQQP MELENKAAVPLGGFLCNVAD KSGAMEMAGLCPAACMQTPR MKLAVQFTNRNQYCYGSRDL LGLHNMKRRQLARLG YRVV



Western blot analysis of various lysates using FASTKD5 Rabbit polyclonal antibody (STJ119787) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000896) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081