

Anti-ZNF768 antibody (140-270) (STJ119807)

STJ119807

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000
	ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79724
Gene Symbol	ZNF768
Uniprot ID	ZN768_HUMAN
Immunogen	
Immunogen Region	140-270
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 140-270 of human ZNF768 (NP_078947.3).
Immunogen Sequence	PGYESESSRYESQNTLKTQ SPEFEAQSSKFQEGAEMLLN PEEKSPLNISVGVHPLDSFT QGFGEQPTGDLPIGPPFEMP TGALLSTPQFEMLQNPLGLT GALRGPGRRGGRARGGGQGPR PNICGICGKSF



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081