

Anti-DPY30 antibody (1-99) (STJ119825)

STJ119825

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84661
Gene Symbol	DPY30
Uniprot ID	DPY30_HUMAN
Immunogen	
Immunogen Region	1-99
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-99 of human DPY30 (NP_115963.1).
Immunogen Sequence	MEPEQMLEGQTQVAENPHSE YGLTDNVERIVENEKINAEK SSKQKVDLQSLPTRYLDQT VVPILLQGLAVLAKERPPNP IEFLASYLLKNKAQFEDRN

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081