

**Anti-HSP90AB1 antibody (1-300) (STJ119988)**  
STJ119988

**GENERAL INFORMATION**


|                          |                        |
|--------------------------|------------------------|
| <b>Product Type</b>      | Primary antibodies     |
| <b>Short Description</b> |                        |
| <b>Applications</b>      | WB/IF/ICC/ELISA        |
| <b>Host/Source</b>       | Rabbit                 |
| <b>Reactivity</b>        | Human/Mouse/Rat/Monkey |

**PRODUCT PROPERTIES**

|                            |                                                                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                                                              |
| <b>Clone ID</b>            |                                                                                                                                                         |
| <b>Concentration</b>       | Lot specific                                                                                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                            |
| <b>Purification</b>        | Affinity purification                                                                                                                                   |
| <b>Dilution</b>            | WB:1:500-1:2000                                                                                                                                         |
| <b>Range</b>               | IF/CC:1:50-1:200<br>ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements. |
| <b>Formulation</b>         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                        |
| <b>Isotype</b>             | IgG                                                                                                                                                     |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                          |

**TARGET INFORMATION**

|                           |                                                                                                                                                                                                                                                                   |
|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 3326                                                                                                                                                                                                                                                              |
| <b>Gene Symbol</b>        | HSP90AB1                                                                                                                                                                                                                                                          |
| <b>Uniprot ID</b>         | HS90B_HUMAN                                                                                                                                                                                                                                                       |
| <b>Immunogen</b>          |                                                                                                                                                                                                                                                                   |
| <b>Immunogen Region</b>   | 1-300                                                                                                                                                                                                                                                             |
| <b>Specificity</b>        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Hsp90 Beta (NP_031381.2).                                                                                                                                            |
| <b>Immunogen Sequence</b> | MPEEVHHGEEVETFAFQAE IAQLMSLIINTFYSNKEIFL RELISNASDALDKIRYESLT DPSKLD SGKELKIDIIPNPQ<br>ERTLTLDVTGIGMTKADLIN NLGTIAKSGTKAFMEALQAG ADISMIGQFGVGFYSAYLVA EKVVVITKHNDDEQYAWESS<br>AGGSFTVRADHGEPIGRGTK VILHLKEDQTEYLEERRVKE VKKKHSQFIGYPITLYLEKE REKEISDDEAEEEEKGEKEE |



Western blot analysis of lysates from wild type (WT) and HSP90AB1 knockout (KO) HeLa cells, using [KO Validated] Hsp90 Beta Rabbit polyclonal antibody (STJ119988) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJSD00856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081