

Anti-ASIC1 antibody (100-150 aa) (STJ13100008)

STJ13100008

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ASIC1 (100-150 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Rat/Mouse/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	41
Gene Symbol	ASIC1
Uniprot ID	ASIC1_HUMAN
Immunogen	A synthetic peptide from aa region 100-150 of the extracellular domain of rat ASIC1 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in mouse and human.
Immunogen Region	100-150 aa
Specificity	Specific for ASIC1.
Immunogen Sequence	