

Anti-Calciressin 1 antibody (STJ13100081)

STJ13100081

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Calciressin 1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A concentration of 10-50 Åµg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	1827
Gene Symbol	RCAN1
Uniprot ID	RCAN1_HUMAN
Immunogen	A synthetic peptide from human Calciressin 1 conjugated to blue carrier protein was used as the antigen. The peptide shares 97% identity with rat and mouse sequence.
Immunogen Region	
Specificity	Specific for Calciressin 1.
Immunogen Sequence	