

Anti-SH3GL2 antibody (C-Term) (STJ13100126)

STJ13100126

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SH3GL2 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	51100
Gene Symbol	SH3GLB1
Uniprot ID	SHLB1_HUMAN
Immunogen	A synthetic peptide from the c-terminal region of human Endophilin-A1 (Endophilin-1) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for SH3GL2.
Immunogen Sequence	