

Anti-HOMER3 antibody (Internal Region) (STJ13100153)

STJ13100153

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HOMER3 (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:2000 is recommended for WB and 1:1000 for IHC-P. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	26558
Gene Symbol	Homer3
Uniprot ID	HOME3_MOUSE
Immunogen	A synthetic peptide from the internal region of rat HOMER3 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	Internal Region
Specificity	Specific for Homer3.
Immunogen Sequence	