

Anti-NDFIP2 antibody (N-Term) (STJ13100191)

STJ13100191

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NDFIP2 (N-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 500 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	Ndfip2
Uniprot ID	NFIP2_MOUSE
Immunogen	A synthetic peptide from the n-terminal of mouse NEDD4 WW domain-binding protein 5A (Ndfip2, NFIP2) conjugated to an immunogenic carrier protein was used as the immunogen.
Immunogen Region	N-Term
Specificity	Specific for Ndfip2.
Immunogen Sequence	