

## Anti-P2RY2 antibody (STJ13100227)

STJ13100227

### GENERAL INFORMATION

|                          |                                                                                                                           |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                        |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-P2RY2 is suitable for use in Immunohistochemistry and Western Blot research applications. |
| <b>Applications</b>      | IHC/WB                                                                                                                    |
| <b>Host/Source</b>       | Rabbit                                                                                                                    |
| <b>Reactivity</b>        | Mouse/Rat                                                                                                                 |

### PRODUCT PROPERTIES

|                       |                                                                                                                                                |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>      | Polyclonal                                                                                                                                     |
| <b>Clone ID</b>       |                                                                                                                                                |
| <b>Concentration</b>  |                                                                                                                                                |
| <b>Conjugation</b>    | Unconjugated                                                                                                                                   |
| <b>Purification</b>   | IgG purified                                                                                                                                   |
| <b>Dilution Range</b> | IHC, WB. A dilution of 1:1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| <b>Formulation</b>    | Lyophilised                                                                                                                                    |
| <b>Isotype</b>        | IgG                                                                                                                                            |
| <b>Storage</b>        | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.          |
| <b>Instruction</b>    | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.                                    |

### TARGET INFORMATION

|                           |                                                                                                  |
|---------------------------|--------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">18442</a>                                                                            |
| <b>Gene Symbol</b>        | <a href="#">P2ry2</a>                                                                            |
| <b>Uniprot ID</b>         | <a href="#">P2RY2_MOUSE</a>                                                                      |
| <b>Immunogen</b>          | A synthetic peptide from mouse P2RY2 conjugated to blue carrier protein was used as the antigen. |
| <b>Immunogen Region</b>   |                                                                                                  |
| <b>Specificity</b>        | Specific for P2RY2.                                                                              |
| <b>Immunogen Sequence</b> |                                                                                                  |