

Anti-POC1A antibody (300-350 aa) (STJ13100239)

STJ13100239

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-POC1A (300-350 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	70235
Gene Symbol	Poc1a
Uniprot ID	POC1A_MOUSE
Immunogen	A synthetic peptide from aa region 300-350 of mouse POC1A conjugated to blue carrier protein was used as the antigen.
Immunogen Region	300-350 aa
Specificity	Specific for POC1A.
Immunogen Sequence	