

Anti-RGS19 antibody (STJ13100295)

STJ13100295

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RGS19 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. Use at a dilution of 1:1000 to 1: 3000. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	10287
Gene Symbol	RGS19
Uniprot ID	RGS19_HUMAN
Immunogen	A synthetic peptide from human RGS19 conjugated to blue carrier protein has been used as the antigen. The peptide is homologous in rat and mouse.
Immunogen Region	
Specificity	Specific for RGS19.
Immunogen Sequence	