

Anti-RAC1+RAC2 antibody (100-150 aa) (STJ13100299)

STJ13100299

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAC1+RAC2 (100-150 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:1000 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	5879
Gene Symbol	RAC1
Uniprot ID	RAC1_HUMAN
Immunogen	A synthetic peptide from aa region 100-150 of human RAC1 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in rat, mouse, guinea pig, bovine, xenopus, dog, zebrafish and chicken.
Immunogen Region	100-150 aa
Specificity	Specific for RAC1 and RAC2.
Immunogen Sequence	